

AG119 Media Replacement Kit for the AG097 Odor Filter

-The first step, remove the AG097 Odor Filter from the suction unit. The next steps are best performed outside, or in an area where carbon dust will not affect other people in the room.

-The top must be removed, by first, removing the screws that hold it to the box.

-The consumed odor material can be emptied into a strong plastic bag or box, sealed, and disposed of properly.

-Next, the first of the two pads should be put into the bottom of the AG097, and then the new odor media poured into the AG097, on top of the pad.

-When full, tap the box gently on a solid surface, to pack the media in well, and then place the second pad directly on the top the new media. -Finally, replace the AG097 cover, and secure it with the screws.

Repl. Media AG119



Remove screws



Pads



Media



Empty AG097



Empty AG097

